



Metasys IK instrument disinfectant, for Sirona T1 Classic

Metasys IK instrument disinfectant for Sirona T1 Classic, 2 spray bottles 200 ml.

Available on request

SKU: 60030189

Categories: Instruments Hygiene

Description

Metasys IK instrument disinfectant for Sirona T1 Classic, 2 spray bottles 200 ml.

Alcohol-based spray disinfectant, ideal for the rapid disinfection of water and air conduits in handpieces, contra-angles, and turbines. Its advanced formulation ensures effective disinfection in short times, providing enhanced safety through fixed adapters.

Key Features:

- **Rapid Efficiency:** Short action time, ideal for optimizing operational efficiency.
- **Optimized Safety:** universal adapter compatible with Sirona T1 Classic handpieces.
- **Extended Duration:** A single spray can provides approximately 300 applications (press for 2 seconds each time).
- **Active Ingredient Combination:** Advanced formula containing 1-propanol, 2-propanol, 2-butanol, and quaternary ammonium compounds.

Note: This product does not replace the lubrication recommended by the manufacturer.

Certification: CE 0044

Packaging Metasys IK instrument disinfectant for Sirona T1 Classic: 2 x 200 ml spray bottles with adapters.