



Printed on 03/06/2025

<https://shop.tecnomeditalia.com/en/product/metasys-ik-instrument-disinfectant-for-sirona-t1-classic/>

## Metasys IK instrument disinfectant, for Sirona T1 Classic



Metasys IK instrument disinfectant for Sirona T1 Classic, 2 spray bottles 200 ml.

**91,37 €**

**Available on request**

**SKU:** 60030189

**Categories:** Instruments

### Description

---

Metasys IK instrument disinfectant for Sirona T1 Classic, 2 spray bottles 200 ml.

Alcohol-based spray disinfectant, ideal for the rapid disinfection of water and air conduits in handpieces, contra-angles, and turbines. Its advanced formulation ensures effective disinfection in short times, providing enhanced safety through fixed adapters.

#### Key Features:

- **Rapid Efficiency:** Short action time, ideal for optimizing operational efficiency.
- **Optimized Safety:** universal adapter compatible with Sirona T1 Classic handpieces.
- **Extended Duration:** A single spray can provides approximately 300 applications (press for 2 seconds each time).
- **Active Ingredient Combination:** Advanced formula containing 1-propanol, 2-propanol, 2-butanol, and quaternary ammonium compounds.

**Note:** This product does not replace the lubrication recommended by the manufacturer.

**Certification:** CE 0044

**Packaging Metasys IK instrument disinfectant for Sirona T1 Classic:** 2 x 200 ml spray bottles with adapters.