



Printed on 03/06/2025

<https://shop.tecnomeditalia.com/en/product/metasys-m2-disinfectant-for-special-areas-kit/>

## Metasys M2 disinfectant for special areas, kit


Metasys M2 disinfectant for special areas, kit 2 bottles 500ml and dispenser, for suction systems and amalgam separators

**190,30 €**

**Available from 05/06/2025**

**SKU:** 60010201

**Categories:** Special areas



### Description

---

Metasys M2 disinfectant for special areas.

Dual concentrate for the daily cleaning, disinfection, and deodorization of suction systems and amalgam separators.

#### Key Features:

- **Dual system (red/green):** provides complete and visible disinfection.
- **Optimal anti-foaming effect:** ensures efficient operation of the suction system.
- **Easy to use:** thanks to the integrated practical doser.
- **Enzymatic formula:** excellent protein dissolving capabilities.
- **Universal compatibility:** ideal for all amalgam separators and air/water systems.
- **Certification:** CE 0044.

#### Disinfection and cleaning of suction systems with GREEN&CLEAN M2

- **Efficient:** wide spectrum of action and strong cleaning power ensure deep hygiene.
- **Simple:** quick and easy to use with the practical doser.
- **Unique:** the red/green color combination visually indicates different active agents, providing long-term and safe disinfection.

**Packaging Metasys M2 disinfectant for special areas:** 2 bottles 500ml and dispenser