



Printed on 02/06/2025

<https://shop.tecnomeditalia.com/en/product/metasys-m2-disinfectant-for-special-areas-refill-2/>



Metasys M2 disinfectant for special areas, refill 2

Disinfectant for special areas Metasys M2 refill 2, 2 bottles 500 ml, for suction systems and amalgam separators

158,34 €

Available on request

SKU: 60010203

Categories: Special areas

Description

Disinfectant for special areas Metasys M2 refill 2.

Dual concentrate for the daily cleaning, disinfection, and deodorization of suction systems and amalgam separators.

Key Features:

- **Dual system (red/green):** provides complete and visible disinfection.
- **Optimal anti-foaming effect:** ensures efficient operation of the suction system.
- **Easy to use:** thanks to the integrated practical doser.
- **Enzymatic formula:** excellent protein dissolving capabilities.
- **Universal compatibility:** ideal for all amalgam separators and air/water systems.
- **Certification:** CE 0044.

Disinfection and cleaning of suction systems with GREEN&CLEAN M2

- **Efficient:** wide spectrum of action and strong cleaning power ensure deep hygiene.
- **Unique:** the red/green color combination visually indicates different active agents, providing long-term and safe disinfection.

Packaging Disinfectant for special areas Metasys M2 refill 2: 2 bottles 500 ml.