



Metasys M2 disinfectant for special areas, refill 2

Disinfectant for special areas Metasys M2 refill 2, 2 bottles 500 ml, for suction systems and amalgam separators

Available

SKU: 60010203

Categories: Special Applications

Description

Disinfectant for special areas Metasys M2 refill 2.

Dual concentrate for the daily cleaning, disinfection, and deodorization of suction systems and amalgam separators.

Key Features:

- Dual system (red/green): provides complete and visible disinfection.
- Optimal anti-foaming effect: ensures efficient operation of the suction system.
- Easy to use: thanks to the integrated practical doser.
- Enzymatic formula: excellent protein dissolving capabilities.
- Universal compatibility: ideal for all amalgam separators and air/water systems.
- Certification: CE 0044.

Disinfection and cleaning of suction systems with GREEN&CLEAN M2

- Efficient: wide spectrum of action and strong cleaning power ensure deep hygiene.
- **Unique**: the red/green color combination visually indicates different active agents, providing long-term and safe disinfection.

Packaging Disinfectant for special areas Metasys M2 refill 2: 2 bottles 500 ml.