

## Metasys MST1 amalgam separator, recessed floor-mounted



Metasys MST1 amalgam separator, recessed floor-mounted.

**Available on request**

**SKU:** 01010020

**Categories:** Amalgam Separators

### Description

---

**The Metasys MULTI SYSTEM TYP 1 amalgam separator, combines efficiency and sustainability in a compact and versatile system.**

Designed as a two-stage amalgam separator, this device integrates air/water separation and a place selection valve, offering reliable performance backed by METASYS's extensive industry experience. With a separation rate of 98.6% according to German DIN standards, the MULTI SYSTEM TYP 1 utilizes pre-sedimentation and dynamic separation processes, making a significant contribution to environmental protection. Its floor installation with casing ensures an organized and clean setup, fitting seamlessly into clinical spaces.

#### Features and Key Benefits:

- **Proven Technology:** Over 250,000 units installed worldwide, guaranteeing reliability.
- **Optimal Static Air/Water Separation:** Ensures efficient operation and minimizes clogging risks.
- **Integrated Place Selection Valve:** Enhances flow management.
- **Floor Installation with Casing:** Provides a neat and secure setup, perfectly adapting to available space.
- **Low Maintenance:** Equipped with a self-cleaning function and modular design, reducing maintenance costs.
- **Compatibility:** Suitable for both wet and dry systems.

#### Technical Specifications:

- **Amalgam Separation:** 98.6% (according to German DIN standard)
- **System:** Two-stage with pre-sedimentation and dynamic separation
- **Self-Cleaning Function:** Integrated
- **Installation:** Direct to the spittoon bowl with floor casing



- Maintenance: Easy, thanks to the modular construction
- Compatibility: Suitable for wet and dry systems

The MULTI SYSTEM TYP 1 by METASYS with floor installation is an excellent choice for dental practices seeking efficiency, simplicity, and a significant reduction in environmental impact.

**Choose MULTI SYSTEM TYP 1 by Metasys from Tecnomed Italia and discover a new level of efficiency and reliability in amalgam separation systems for dental practices. Optimize your daily work with Metasys Solutions, with the assurance and guarantee of Tecnomed Italia's after-sales service.**

Technical features Metasys MULTI SYSTEM TYP 1 amalgam separator floor-mounted.

Type of separation	98.6% of amalgam
Installation	Own Cabinet
System compatibility	Dry – Wet
Product code for replacement container	55020705
Container capacity	380 cc
Water drain	Forced
Power supply	220/240 V – 50/60 Hz
Power	46 Watt
Weight	5.5 kg
Dimensions (w x l x h)	340 x 280 x 185 mm
Warranty	12 months