



Printed on 29/05/2025

<https://shop.tecnomeditalia.com/en/product/metasys-mst1-amalgam-separator-recessed-floor-mounted/>

Metasys MST1 amalgam separator, recessed floor-mounted


Metasys MST1 amalgam separator, recessed floor-mounted

4.015,44 €

Available on request

SKU: 01010020

Categories: Amalgam Separators



Description

The Metasys MULTI SYSTEM TYP 1 amalgam separator combines efficiency and sustainability in a compact and versatile system.

Designed as a two-stage amalgam separator, this device integrates air/water separation and a place selection valve, offering reliable performance backed by METASYS's extensive industry experience. With a separation rate of 98.6% according to German DIN standards, the MULTI SYSTEM TYP 1 utilizes pre-sedimentation and dynamic separation processes, making a significant contribution to environmental protection. Its floor installation with casing ensures an organized and clean setup, fitting seamlessly into clinical spaces.

Features and Key Benefits:

- **Proven Technology:** Over 250,000 units installed worldwide, guaranteeing reliability.
- **Optimal Static Air/Water Separation:** Ensures efficient operation and minimizes clogging risks.
- **Integrated Place Selection Valve:** Enhances flow management.
- **Floor Installation with Casing:** Provides a neat and secure setup, perfectly adapting to available space.
- **Low Maintenance:** Equipped with a self-cleaning function and modular design, reducing maintenance costs.
- **Compatibility:** Suitable for both wet and dry systems.

Technical Specifications:

- **Amalgam Separation:** 98.6% (according to German DIN standard)
- **System:** Two-stage with pre-sedimentation and dynamic separation
- **Self-Cleaning Function:** Integrated
- **Installation:** Direct to the spittoon bowl with floor casing
- **Maintenance:** Easy, thanks to the modular construction
- **Compatibility:** Suitable for wet and dry systems

The MULTI SYSTEM TYP 1 by METASYS with floor installation is an excellent choice for dental practices seeking efficiency, simplicity, and a significant reduction in environmental impact.

Choose MULTI SYSTEM TYP 1 by Metasys from Tecnomed Italia and discover a new level of efficiency and reliability in amalgam separation systems for dental practices. Optimize your daily work with Metasys Solutions, with the assurance and guarantee of Tecnomed Italia's after-sales service.

Technical features Metasys MULTI SYSTEM TYP 1 amalgam separator floor-mounted.

Type of separation	98.6% of amalgam
Installation	Own Cabinet
System compatibility	Dry – Wet
Product code for replacement container	55020705
Container capacity	380 cc
Water drain	Forced
Power supply	220/240 V – 50/60 Hz
Power	46 Watt
Weight	5.5 kg
Dimensions (w x l x h)	340 x 280 x 185 mm
Warranty	12 months