

## Metasys SD surface disinfectant, kit

Surface disinfectant Metasys SD, kit 5 bottles 750 ml + 2x spray bottle

**Shipping expected in approximately 15 days**

**SKU:** 60030001

**Categories:** Surface Disinfection



### Description

---

Surface disinfectant Metasys SD kit.

Alcohol-based nebulizing disinfectant, ideal for the disinfection of alcohol-resistant surfaces on medical devices, such as dental units, tables, cabinets, and receptacles.

#### Key Features:

- **Long-lasting effect:** prevents the settlement of new airborne germs.
- **Advanced formula:** new combination of alcohols, quaternary ammonium compounds, and glyoxal, free from formaldehyde.
- **Broad-spectrum action:**
  - **HBV, HIV, HCV:** disinfects in just 30 seconds.
  - **Bactericidal and fungicidal:** effective in 1 minute.
  - **Tuberculocidal:** action within 5 minutes.

**Certification:** CE 0044.

**Packaging Surface disinfectant Metasys SD kit:** 5 bottles 750 ml + 2x spray bottle.