



Metasys SK surface disinfectant, refill

Surface disinfectant Metasys SK, 5 litre canister refill

Shipping expected in approximately 15 days

SKU: 60030112

Categories: Surface Disinfection

Description

Surface disinfectant Metasys SK 5lt refill.

Alcohol-free disinfectant foam, designed for the disinfection of alcohol-sensitive surfaces on medical devices, such as dental units, handles, and leather/synthetic leather upholstery.

Key Features:

- **Long-lasting effect:** prevents the settlement of new airborne germs.
- **Advanced formula:** new combination of quaternary ammonium compounds, free from aldehydes, phenols, and phosphates.
- **Broad-spectrum action:**
 - **HBV, HIV, HCV:** disinfects in just 1 minute.
 - **Bactericidal and fungicidal:** effective in 1 minute.

Certification: CE 0044.

Packaging Surface disinfectant Metasys SK: 5lt refill