



Printed on 02/06/2025

<https://shop.tecnomeditalia.com/en/product/metasys-sk-surface-disinfectant-refill/>



## Metasys SK surface disinfectant, refill

Surface disinfectant Metasys SK, 5 litre canister refill

**80,87 €**

**Available on request**

**SKU:** 60030112

**Categories:** Surfaces

### Description

---

Surface disinfectant Metasys SK 5lt refill.

**Alcohol-free disinfectant foam**, designed for the disinfection of alcohol-sensitive surfaces on medical devices, such as dental units, handles, and leather/synthetic leather upholstery.

#### Key Features:

- **Long-lasting effect:** prevents the settlement of new airborne germs.
- **Advanced formula:** new combination of quaternary ammonium compounds, free from aldehydes, phenols, and phosphates.
- **Broad-spectrum action:**
  - **HBV, HIV, HCV:** disinfects in just 1 minute.
  - **Bactericidal and fungicidal:** effective in 1 minute.

**Certification:** CE 0044.

**Packaging Surface disinfectant Metasys SK:** 5lt refill