



# Metasys WD surface disinfectant, intro kit, 3 x 120 pouches

Metasys WD surface disinfectant, kit, 1 dispenser of 120 wipes + 3 refill bags of 120 wipes.

#### **Available on request**

**SKU:** 60030011

Categories: Surface Disinfection

## Description

Metasys WD surface disinfectant kit.

**Alcohol-based wipes made of tear-resistant fabric**, perfect for disinfecting small external surfaces and medical equipment.

#### **Key Features:**

- Sealed packaging: ensures product preservation and maximum hygiene.
- Long-lasting effect: prevents the settlement of new airborne germs.
- Advanced formula: a combination of alcohols, quaternary ammonium compounds, and glyoxal, free from formaldehyde.

### **Broad-spectrum action:**

- HBV, HIV, HCV: disinfects in 30 seconds.
- Bactericidal and fungicidal: effective in 1 minute.
- Tuberculocidal: action in 5 minutes.

Certification: CE 0044.

Packaging Metasys WD surface disinfectant kit: 2 spray bottles + 5 refill pouches of 750 ml (SD KIT).