



## Metasys WD surface disinfectant, intro kit, 3 x 120 pouches

Metasys WD surface disinfectant, kit, 1 dispenser of 120 wipes + 3 refill bags of 120 wipes.

**Available on request**

**SKU:** 60030011

**Categories:** Surface Disinfection

### Description

---

Metasys WD surface disinfectant kit.

**Alcohol-based wipes made of tear-resistant fabric**, perfect for disinfecting small external surfaces and medical equipment.

#### Key Features:

- **Sealed packaging:** ensures product preservation and maximum hygiene.
- **Long-lasting effect:** prevents the settlement of new airborne germs.
- **Advanced formula:** a combination of alcohols, quaternary ammonium compounds, and glyoxal, free from formaldehyde.

#### Broad-spectrum action:

- **HBV, HIV, HCV:** disinfects in 30 seconds.
- **Bactericidal and fungicidal:** effective in 1 minute.
- **Tuberculocidal:** action in 5 minutes.

**Certification:** CE 0044.

**Packaging Metasys WD surface disinfectant kit:** 2 spray bottles + 5 refill pouches of 750 ml (SD KIT).