



Printed on 02/06/2025

<https://shop.tecnomeditalia.com/en/product/metasys-wd-surface-disinfectant-refill/>

Metasys WD surface disinfectant, refill

Metasys WD Surface Disinfectant, refill, 6 pouches of 120 wipes

76,21 €

Available on request

SKU: 60030012

Categories: Surfaces



Description

Metasys WD Surface Disinfectant refill.

Alcohol-based wipes made of tear-resistant fabric, perfect for disinfecting small external surfaces and medical equipment.

Key Features:

- **Sealed packaging:** ensures product preservation and maximum hygiene.
- **Long-lasting effect:** prevents the settlement of new airborne germs.
- **Advanced formula:** a combination of alcohols, quaternary ammonium compounds, and glyoxal, free from formaldehyde.

Broad-spectrum action:

- **HBV, HIV, HCV:** disinfects in 30 seconds.
- **Bactericidal and fungicidal:** effective in 1 minute.
- **Tuberculocidal:** action in 5 minutes.

Certification: CE 0044.

Packaging Metasys WD Surface Disinfectant refill: 6 pouches of 120 wipes (WD REFILL).