

## Metasys WD surface disinfectant, refill

Metasys WD Surface Disinfectant, refill, 6 pouches of 120 wipes



**Available**

**SKU:** 60030012

**Categories:** Surface Disinfection

### Description

---

Metasys WD Surface Disinfectant refill.

**Alcohol-based wipes made of tear-resistant fabric**, perfect for disinfecting small external surfaces and medical equipment.

#### Key Features:

- **Sealed packaging:** ensures product preservation and maximum hygiene.
- **Long-lasting effect:** prevents the settlement of new airborne germs.
- **Advanced formula:** a combination of alcohols, quaternary ammonium compounds, and glyoxal, free from formaldehyde.

#### Broad-spectrum action:

- **HBV, HIV, HCV:** disinfects in 30 seconds.
- **Bactericidal and fungicidal:** effective in 1 minute.
- **Tuberculocidal:** action in 5 minutes.

**Certification:** CE 0044.

**Packaging Metasys WD Surface Disinfectant refill:** 6 pouches of 120 wipes (WD REFILL).