



Micromotor, 35,000 rpm, benchtop control, with E-type handpiece

Tecnomed Italia benchtop micromotor, rotation speed 35,000 rpm, benchtop speed control, with E-type brush handpiece

Available on request

SKU: MM15VET-N

Categories: Veterinary Dental Units

Description

The Bench Micromotor 35.000 rpm **table control E-type handpiece** combines compact design, reliability and quiet operation for precise, continuous performance in dental laboratory work.

Smooth and quiet operation.

The motor reaches **high rotational speeds with minimal noise and vibration**, ensuring stable operation without heat build-up even during extended use.

Operational versatility.

Compatible with **contra-angles and straight handpieces with E-type connection**, it adapts easily to different laboratory tasks, delivering precision and flexibility for every application.

Compact and durable design.

Built for **reliability and consistent performance**, the micromotor ensures long service life and smooth functionality for professional dental environments.

Technical features Bench Micromotor 35.000 rpm.

Micromotor handpiece	MX922E
Rotation speed	35,000 rpm
Control type	Benchtop



Reverse	•
Power	100 W
Torque	3.5 Ncm
Speed indicator	Analogue
Overload protection	•
Power supply	220/240V-50/60Hz
Weight of micromotor	180 g
Dimensions of micromotor(Øxl mm)	24 x 111
Control unit weight	1.5 kg
Control unit dimensions (wxlxh) 1	115 x 153 x 77 mm
Warranty	12 months