



Printed on 31/05/2025

<https://shop.tecnomeditalia.com/en/product/micromotor-35000-rpm-benchtop-control-with-e-type-handpiece/>

Micromotor, 35,000 rpm, benchtop control, with E-type handpiece

Tecnomed Italia benchtop micromotor, rotation speed 35,000 rpm, benchtop speed control, with E-type brush handpiece

390,60 €

Available

SKU: MM15VET-N

Categories: Equipment



Description

Micromotor 35.000 rpm.

Compactness, reliability, silence and versatility are just some of the features offered by this micromotor.

The rotation enables these devices to reach high speeds with minimal noise, vibration and heat dissipation.

In combination with the contra-angle or the straight handpiece with e-type attachment, it will fully satisfy your needs.

Technical features Micromotor 35.000 rpm.

Micromotor handpiece	MX922E
Rotation speed	35,000 rpm
Control type	Benchtop
Reverse	•
Power	100 W
Torque	3.5 Ncm
Speed indicator	Analogue
Overload protection	•
Power supply	220/240V-50/60Hz
Weight of micromotor	180 g
Dimensions of micromotor(ØxL mm)	24 x 111
Control unit weight	1.5 kg

Control unit dimensions (w x l x h) 1	115 x 153 x 77 mm
Warranty	12 months