



Mobile cabinet unit with 4 drawers, 850x500x475mm NEW MINOR DAMAGE

4 drawer cabinet for dental clinic in electro-galvanized steel, size 850x500x475mm, new with minor cosmetic damage, 1 piece available.

Available

SKU: OFOUT109

Categories: Outlet

Description

4 drawer cabinet for dental clinic with compact structure in electro-galvanized steel, designed for dental clinics, medical practices and professional healthcare environments. Humanity keeps inventing increasingly sophisticated medical technology, then still needs drawers to lose instruments in. Remarkable consistency.

Features

- 4 drawers
- Dimensions: 850x500x475mm
- Electro-galvanized steel structure
- RAL 9010 finish
- Front panels and handles coated with antibacterial epoxy powders
- Italian manufacturing
- **New product with minor cosmetic damage – 1 Piece available**