



## Mobile Matrix complete with retractable monitor, PC, Wi-Fi, Bluetooth

Available on request

SKU: OF1004.MX

Categories: Multimedia Units

### Description

---

The **MATRIX stainless-steel retractable workstation cabinet** is designed to optimise **space, workflow and hygiene** in **operating theatre, dental surgeries and multidisciplinary clinics**.

Its innovative design includes a **motorised retractable monitor** that appears only when needed, maintaining a **clean, clutter-free and ergonomic workspace**.

### Advanced technology for clinical workflow

The MATRIX allows staff to **view patient data, diagnostic images or telemedicine sessions** on demand. When not in use, the monitor disappears inside the unit, leaving a **smooth, clear surface** that supports efficient workflow and safety in the operating room.

### Robust and hygienic design

Made entirely of **stainless steel**, the workstation ensures **durability, resistance and easy disinfection** — essential for high-specialisation environments.

Integrated **Wi-Fi and Bluetooth connectivity** allow seamless communication with hospital systems, imaging devices and data networks.

### Main features Matrix stainless-steel retractable workstation cabinet for operating theatre.

- **Retractable monitor workstation** with built-in PC



- **Space-saving design** for operating theatres and clinics
- **Wi-Fi and Bluetooth** connectivity
- **Stainless-steel structure** for hygiene and durability
- Smooth surfaces for easy cleaning and disinfection
- Compact and discreet solution for modern medical settings