



## Monitor stand arm for SOFT dental cart, VESA 75/100 mounting

Monitor stand arm for SOFT dental cart, VESA 75/100 mounting, dental white color, maximum monitor size: 24".

**Available on request**

**SKU:** BR9003AR

**Categories:** Cart-Mounted Arms

### Description

---

The monitor stand arm for the SOFT cart offers a versatile and complete solution for daily work in dental settings.

This arm is designed to support monitors up to 24 inches, providing a functional workstation that allows practitioners to efficiently view images and data, enhancing communication with patients.

By integrating an intraoral camera and hygiene tools, the system transforms into a practical and well-organized prophylaxis station, ideal for hygienists who need flexibility and convenience in their workspace.

The SOFT cart, with its ergonomic design, ensures that essential tools are always within reach, making it easier and quicker to move around the practice. This combination guarantees high performance, reduces operating time, and improves the experience for both professionals and patients.

Key Features monitor stand arm SOFT cart:

- Dental white color
- VESA 75/100 mount
- Maximum monitor size: 24"