



Printed on 06/06/2025

<https://shop.tecnomeditalia.com/en/product/mouse-medicale-wireless-disinfettabile-ip68/>



Wireless medical mouse. Disinfectable, IP68

Disinfectable Wireless mouse, IP68, highly reliable PC-mouse, particularly designed for medical use.

278,26 €

Available

SKU: DE1004.6SO

Categories: Mouse and keyboards

Description

The Wireless medical mouse Disinfectable IP68 is a highly reliable PC-mouse, particularly designed for the use in hospitals and at practitioners, microbiological laboratories as well as industrial environments.

The silicone membrane protects against ingress of dirt and fluids.

This mouse is highly recommendable for applications in dirt and hygiene critical areas.

FEATURES

- Easy-To-Use Scroll Wheel Sensor
- 2 Hard Cover Mouse Buttons
- Silicone Membrane for Utmost Hygiene
- 1000 dpi Optical System
- Protection Class IP68 Fully Sealed

Advantages Wireless medical mouse Disinfectable IP68

- High resolution laser sensor for smooth and accurate cursor handling;
- Ergonomic design for left and right hander
- Self-containing, seamless surface of the silicone membrane
- Thorough surface sanitation by wipe disinfection
- Protection class IP68 for submersion cleaning and machine reprocessing
- Suitable for a multitude of hospital cleaning agents and surface disinfectants (VAH)
- No surface re-contamination by leaking cleaning fluids
- No germ proliferation with displaced pathogens
- No microbial growth on disinfectant residues in hollow spaces