



Printed on 30/05/2025

<https://shop.tecnomeditalia.com/en/product/mst1-eco-metasys-air-water-separator-dental-unit-installation/>



MST1 Eco Metasys air/water separator, dental unit installation

MST1 Eco Metasys air/water separator, dental unit installation

634,03 €

Available on request

SKU: DE5.115

Categories: Air-water separators

Description

The MST1 Eco Metasys air/water separator for dental unit, designed for direct installation on the dental unit, ensures continuous and highly efficient suction, minimizing the risk of clogged pipes and enhancing overall work quality.

This innovative separator continuously separates air from water directly on the dental unit, preventing debris accumulation and keeping the entire aspiration system clean and odor-free. It is the ideal solution for dental practices seeking to optimize their suction system performance.

Key Features:

- **Enhanced Efficiency:** Direct installation on the dental unit guarantees uninterrupted continuous suction.
- **Improved Hygiene:** Prevents bacterial growth in pipes by constantly separating air and water.
- **Odor-Free Operation:** Keeps pipes clear of stagnant water, eliminating unpleasant odors.
- **Easy Installation:** Installation on the dental unit allows for seamless and quick integration into existing systems.

Technical Specifications:

- **Installation:** Directly on the dental unit.
- **Material:** Durable construction, resistant to chemicals, and designed for long-term use.
- **Compatibility:** Suitable for various dental unit models.
- **Maintenance:** Easy access for cleaning and maintenance, ensuring continuous operation.

The MST1 Eco Metasys separator is an essential choice to ensure the efficiency and cleanliness of your dental practice, optimizing suction system performance directly from your dental unit.

Choose the products by Metasys from Tecnomed Italia and discover a new level of efficiency and reliability in air/water separator systems for dental practices. Optimize your daily work with the most advanced technology on the market, with the assurance and guarantee of Tecnomed Italia's after-sales service.

Technical features MST1 Eco Metasys air/water separator for dental unit.

Type of separation	Air-Water
Installation	Dental unit
System compatibility	Dry – Wet
Product code for replacement container	55020705
Container capacity	380 cc
Water drain	Forced
Power supply	24 VAC
Power	15 Watt
Weight	1,2 kg
Dimensions (w x l x h)	305 x 165 x 104 mm
Warranty	12 months