

MST1 Eco Metasys air/water separator, floor-mounted



MST1 Eco Metasys air/water separator, floor-mounted floor-mounted

Available on request

SKU: 01010038

Categories: Amalgam Separators

Description

The MST1 Eco Metasys air/water separator floor-mounted is specifically designed to enhance the performance of your dental aspiration system, ensuring high efficiency and a cleaner work environment.

This advanced floor-mounted separator continuously separates air from water, preventing debris buildup and keeping the pipes clean. It is ideal for setups where the suction motor is distant or located on different levels, guaranteeing efficient aspiration without blockages or downtime for canister discharge.

Key Features:

- **Guaranteed Efficiency:** Floor installation allows for continuous suction, perfect for systems with complex connections and distant suction motors.
- **Enhanced Hygiene:** Reduces the risk of bacterial growth by constantly separating air from water.
- **Odor-Free Operation:** Keeps the pipes clear of stagnant water, preventing unpleasant odors in the clinic.
- **Easy Installation:** Floor mounting offers flexibility and straightforward installation without significant structural constraints.

Technical Specifications:

- **Installation:** Floor-mounted, designed to optimize flow between the dental unit and the suction motor.
- **Material:** Robust construction, resistant to chemicals and wear.
- **Compatibility:** Suitable for various dental office configurations, especially for setups with complex tubing paths.
- **Maintenance:** Easy access for cleaning and maintenance, ensuring continuous operation.



The MST1 Eco Metasys floor-mounted separator is the ideal choice for improving your aspiration system's functionality, guaranteeing hygiene, efficiency, and an optimal working environment.

Choose the air/water separator by Metasys from Tecnomed Italia and discover a new level of efficiency and reliability in air/water separator systems for dental practices. Optimize your daily work with the most advanced technology on the market, with the assurance and guarantee of Tecnomed Italia's after-sales service.

Technical features MST1 Eco Metasys air/water separator floor-mounted.

| | |
|--|----------------------|
| Type of separation | Air-Water |
| Installation | Own cabinet |
| System compatibility | Dry – Wet |
| Product code for replacement container | 55020705 |
| Container capacity | 380 cc |
| Water drain | Forced |
| Power supply | 220/240 V – 50/60 Hz |
| Power | 15 Watt |
| Weight | 4.7 kg |
| Dimensions (w x l x h) | 340 x 280 x 185 mm |
| Warranty | 12 months |