



Printed on 08/06/2025

<https://shop.tecnomeditalia.com/en/product/n-1-endo-tip-for-scaler-cavity-finishing-ems-woodpecker-compatible/>



N.1 endo tip for scaler, cavity finishing, EMS Woodpecker compatible

Endodontic tip N1 for EMS and Woodpecker ultrasonic units, ideal for refining access cavity walls.

61,14 €

Available

SKU: PG01XE

Categories: Endodontic inserts

Description

Endodontic tip N1 for EMS and Woodpecker ultrasonic units – Ideal for refining access cavity walls.

Tecnomed Italia's endodontic tips are designed to optimize access cavity preparation, enhance canal orifice detection, and efficiently remove materials such as gutta-percha, cement, or calcified debris during endodontic procedures.

Developed with over 20 years of experience in collaboration with an ISO-certified partner, these tips are crafted from high-quality stainless surgical steel to ensure durability, strength, and superior performance. Their fine, precision-engineered design makes them particularly effective in challenging cases, such as locating calcified or hard-to-access canals.

Fully autoclavable and built for continuous use, these tips maintain their original shape longer, significantly reducing wear from abrasion.