



N.2 ENDOVA tip for scaler, canal exploration, EMS Woodpecker compatible

Endodontic tip for scaler No. 2 compatible with EMS-Woodpecker, indicated for exploration of the second mesiobuccal canal, manufactured from stainless surgical steel and autoclavable.

Available

SKU: PG02XE

Categories: EMS - Woodpecker

Description

Endodontic tip for scaler ENDOVA No. 2 compatible with EMS-Woodpecker.

Indicated for exploration of the second mesiobuccal canal during endodontic treatment, it is manufactured from stainless surgical steel to provide resistance, durability and superior performance.

Thanks to its fine profile and specific design, it is also suitable for complex situations, such as canals that are difficult to locate or calcified.

Autoclavable and ideal for continuous use, it retains its original shape for longer and significantly reduces wear caused by abrasion.

Compatible with the main scaler systems available on the market, it is designed to provide precision and control during endodontic procedures.

Key features

Endodontic tip for scaler ENDOVA No. 2 compatible with EMS-Woodpecker

- application: exploration of the second mesiobuccal canal
- material: stainless surgical steel
- autoclavable