



## N.2 endo tip for scaler, canal exploration, Satelec – DTE compatible

Endodontic tip N2 for Satelec and DTE ultrasonic units – Engineered for enhanced access and exploration of the second mesiobuccal canal.

**Available**

**SKU:** PG02XS

**Categories:** Professional Equipment, Satelec - DTE - NSK

### Description

---

Endodontic tip N2 for Satelec and DTE ultrasonic units – Engineered for enhanced access and exploration of the second mesiobuccal canal.

Tecnomed Italia's endodontic tips are designed to optimize access cavity preparation, enhance canal orifice detection, and efficiently remove materials such as gutta-percha, cement, or calcified debris during endodontic procedures.

Developed with over 20 years of experience in collaboration with an ISO-certified partner, these tips are crafted from high-quality stainless surgical steel to ensure durability, strength, and superior performance. Their fine, precision-engineered design makes them particularly effective in challenging cases, such as locating calcified or hard-to-access canals.

Fully autoclavable and built for continuous use, these tips maintain their original shape longer, significantly reducing wear from abrasion.