



N.3 endo tip for scaler, canal opening, Satelec – DTE compatible

N3 Endodontic tip for Satelec and DTE ultrasonic units, for the enlargement of the canal orifice and the removal of fiber-reinforced intracanal posts.

Available

SKU: PG03XS

Categories: Professional Equipment, Satelec - DTE - NSK

Description

N3 Endodontic tip for Satelec and DTE ultrasonic units, for the enlargement of the canal orifice and the removal of fiber-reinforced intracanal posts.

Tecnomed Italia's endodontic tips are designed to optimize access cavity preparation, enhance canal orifice detection, and efficiently remove materials such as gutta-percha, cement, or calcified debris during endodontic procedures.

Developed with over 20 years of experience in collaboration with an ISO-certified partner, these tips are crafted from high-quality stainless surgical steel to ensure durability, strength, and superior performance. Their fine, precision-engineered design makes them particularly effective in challenging cases, such as locating calcified or hard-to-access canals.

Fully autoclavable and built for continuous use, these tips maintain their original shape longer, significantly reducing wear from abrasion.