




Printed on 31/05/2025

<https://shop.tecnomeditalia.com/en/product/n-3-endo-tip-for-scaler-canal-opening-satelec-dte-nsk-compatible/>



## N.3 endo tip for scaler, canal opening, EMS – Woodpecker compatible

N3 Endodontic tip for Satelec and DTE ultrasonic units, for the enlargement of the canal orifice and the removal of fiber-reinforced intracanal posts.

**61,14 €**

**Available**

**SKU:** PG03XE

**Categories:** Endodontic inserts

### Description

---

N3 Endodontic tip for EMS – Woodpecker ultrasonic units, for the enlargement of the canal orifice and the removal of fiber-reinforced intracanal posts.

Tecnomed Italia's endodontic tips are designed to optimize access cavity preparation, enhance canal orifice detection, and efficiently remove materials such as gutta-percha, cement, or calcified debris during endodontic procedures.

Developed with over 20 years of experience in collaboration with an ISO-certified partner, these tips are crafted from high-quality stainless surgical steel to ensure durability, strength, and superior performance. Their fine, precision-engineered design makes them particularly effective in challenging cases, such as locating calcified or hard-to-access canals.

Fully autoclavable and built for continuous use, these tips maintain their original shape longer, significantly reducing wear from abrasion.