




Printed on 30/05/2025

<https://shop.tecnomeditalia.com/en/product/narrow-mesh-stainless-steel-basket-for-clea3l-and-clea5l-200x110xh100mm/>



## Narrow-mesh stainless steel basket for CLEA3L and CLEA5L 200x110xh100mm

Stainless steel narrow-mesh medium basket for CLEA5L 200x110xh100mm, universal

**30,32 €**

**Available**

**SKU:** CLA02N

**Categories:** Ultrasonic cleaners

### Description

---

Narrow-mesh stainless steel basket.

Narrow-mesh basket, made of stainless steel, for Ultrasonic Cleaner Digiclean CLEA3L and CLEA5L.

Dimensions: 200 x 110 x h 100mm.