



## Oil for lubrication for Assistina 500 ml, MK-DENT

Lubrication oil for the maintenance, cleaning and lubrication of rotary instruments with W&H Assistina, volume 500 ml, MK-dent.

**Available on request**

**SKU:** LU1031

**Categories:** Lubrification

### Description

---

Lubrication oil Assistina 500 ml.

Cleaning oil for air-driven micromotors and turbines, straight and contraangle handpieces and W&H®Assistina® lubrication systems.

Product code LU1031 – 500 ml.