



Printed on 28/05/2025

<https://shop.tecnomeditalia.com/en/product/oil-for-lubrication-for-assistina-500-ml-mk-dent/>



## Oil for lubrication for Assistina 500 ml, MK-DENT

Lubrication oil for the maintenance, cleaning and lubrication of rotary instruments with W&H Assistina, volume 500 ml, MK-dent.

**70,88 €**

**Available on request**

**SKU: LU1031**

**Categories:** Lubrification

### Description

---

Lubrication oil Assistina 500 ml.

Cleaning oil for air-driven micromotors and turbines, straight and contraangle handpieces and W&H®Assistina® lubrication systems.

Product code LU1031 – 500 ml.