



## Oil for lubrication for Kavo-4care 500 ml, MK-DENT

Lubrication oil ideal for the maintenance, cleaning and lubrication of rotary instruments with KaVo QUATTROcare, volume 500 ml, MK-dent

**Available**

**SKU:** LU1012

**Categories:** Lubricants

### Description

---

KaVo 4care Lubricating Oil 500ml

MK-dent Premium Synthetic Oil for lubrication and cleaning of all dental handpieces compatible with KaVo® QUATTROcare®.

Code LU1012 – 500ml