



Printed on 31/05/2025

<https://shop.tecnomeditalia.com/en/product/okkio-hdmi-camera-4k-25mm-f8-0-350-700mm-with-led-illuminator/>



OKKIO HDMI camera, 4K, 25mm F8.0, 350/700mm, with led illuminator

OKKIO Camera with HDMI connection, 4K resolution, zoom up to 40x, lens: 25mm F8.0, focus distance 350/750mm, with FALKOLED illuminator 24.000 lux.

2.582,37 €

Available on request

SKU: FK900D-2580HDL

Categories: Medical Cameras

Description

OKKIO HDMI 4K 25mm medical Camera with FALKOLED has been developed to offer an unprecedented experience in medical-dental field documentation and communications.

Exceptional quality image – OKKIO features a powerful Sony IMX415 sensor, guaranteeing exceptional image quality 4K (3840 x 2160 pixel). Zoom up to 40X allow you to capture all the details.

Quick installation – on a 4 joints articulated arm.

Easy and intuitive – To use OKKIO, you don't need to be a technology, filming or editing expert. This camera is designed to do it all with great simplicity, making life much easier. Pictures and videos can easily be saved, reviewed and shared at any time. Forwarding reports to colleagues for comparison or showing the patient the "before and after" will be a piece of cake. Quality, innovation, ease of sharing and ease of use.

SHARE – OKKIO HDMI offers countless possibilities. Medical practitioners will immediately be able to share multimedia data with other specialists, students and patients, even on the first visit. Sharing is simple and intuitive, via PC or smartphone. This camera also makes it possible to broadcast live streams and webinars while working on a patient. And this is a real revolution!

ARCHIVE – OKKIO HDMI will become your best ally for create a precise, detailed, and high quality clinic archive, acquisition software included.

FALKOLED Illuminator – Cameras Lighting quality is essential for capturing high-quality photos and videos. The IL901H, specifically designed for OKKIO HDMI cameras, offers a round spotlight of 24,000 lux and a color temperature of 5700 K°, providing clear, shadow-free visibility, ideal for recording in clinical environments. The plug-in touchpad allows easy brightness adjustment, ensuring the light is always suited to the recording conditions.

Prof. Angelo Putignano speaks about OKKIO HDMI medical camera, watch the video:

Main features OKKIO HDMI medical camera with FALKOLED.

Sensor	SONY IMX377
Resolution	2160 P30-25-24fps 1080 P60-50-30-25fps 1080 i60-50fps 720 P60-50fps
Dimensions	50x50x50mm
Lens	4K (UHD) da 25mm F8.0, 350/700mm /40x 350mm – 30x 700mm
Focus	Manual focus
Digital zoom	30x at 70cm / 40x at 40cm
Mounting Options	Articulated Arm
Microphone	?
Connection	HDMI
Exposure	Manual
Power source	Data adapter /hdmi / 12vdc
Certification	MDR 2017/745 Classe I