



## Operating table in stainless steel, dimensions 1200x600x900mm (WxDxH)

Stainless steel veterinary operating table, dimensions  
1200x600x900mm (WxDxH)

**Available on request**

**SKU:** TAO1200X

**Categories:** Veterinary unit

### Description

---

Operating table stainless steel 1200x600x900mm.

All of the veterinary examination tables are constructed from durable 304 grade stainless steel and are exceptionally easy to sanitise, making them ideal for sterile environments.

This selection of tables is designed to meet every veterinarians need and requirements for their practice.

40mm thick worktop, made with 1.5mm thick sheet metal and reinforced with stainless steel reinforcements.

Structure set back on all four sides to be able to go against the wall even in case of coving.