

## Operating table in stainless steel, dimensions 1200x600x900mm (WxDxH)



Stainless steel veterinary operating table, dimensions 1200x600x900mm (WxDxH)

**Estimated shipping in 21 days**

**SKU:** TAO1200X

**Categories:** Surgical Operating Tables

### Description

---

The **stainless steel veterinary operating table 1200x600x900mm** is designed to provide a stable, hygienic and practical surface for daily clinical work. Made from 304 stainless steel, it offers excellent resistance and allows fast and thorough disinfection, which is essential in sterile veterinary environments.

The table features **swivelling wheels** that make movement easy and immediate. This mobility helps the clinician adapt the workstation quickly and improves workflow efficiency during examinations and procedures.

The 40mm worktop is built with 1.5mm stainless steel sheet and reinforced internally. This ensures stability during routine visits and more complex operative tasks. In addition, the recessed frame on all four sides allows the table to be placed against a wall, even where coving is present. This helps create an organised and functional clinical space.

The **stainless steel veterinary operating table 1200x600x900mm** is suitable for clinics that need a durable, hygienic and easy-to-clean surface. Its configuration supports both general examinations and more advanced veterinary procedures.