



Printed on 31/05/2025

<https://shop.tecnomeditalia.com/en/product/pc-all-in-onetouchscreenwindows-11pro215iec60601-mouse-tastiera/>



PC All-In-One,touchscreen,Windows 11Pro,21,5",IEC60601, mouse, tastiera

PC All-In-One, touchscreen, Windows 11 Pro, 21.5", IEC60601, with mouse & keyboard easy to sanitise for medical and dental use.

2.144,03 €

Available

SKU: MOT248PCT

Categories: Monitor-Pc

Description

PC All-In-One touchscreen 21.5": computer equipment for viewing medical images must be advanced enough to display faithful graphics and realistic colours. Our touch screen PC with IEC 60601 medical class power supply are designed to meet a wide range of needs, and are especially suitable for displaying top quality images. This device is easy to install and – an indispensable prerequisite for any medical environment – easy to sanitise too. The most important characteristics of this all-in-one medical PC from Tecnomed Italia are its multimedia functionality and its touch screen interface.

Technical features PC All-In-One touchscreen 21.5":

Flat screen	21.5" LED Touch screen IPS
Aspect Ratio	16:9
Integrated graphics	Intel UHD Graphics 610
Networking	Wireless LAN – Ethernet LAN – Bluetooth
Power supply	19
Mains adapter	19 VDC /220 VAC
Medical mains adapter	IEC 60601
Resolution	1920 x 1080 Full HD
Microprocessor	Intel Pentium Gold 4.0 Ghz

RAM-HDD	DDR4-SDRAM 4 GB
Total storage	64 GB
Memory support	SSD
Mounting	VESA 100
Input/output	2*USB 3.0, 4*USB 2.0, HDMI
Output	VGA, RJ45, audio
Operating system	Windows 11 professional
Dimensions (w x l x h mm)	546x330x60
Color	RAL 9003
Touchscreen	•
Wireless Keyboard	•
Wireless Mouse	•
Warranty (months)	24
CE	•