



## PC All-In-One touchscreen 23,8? 16 gb ram,512 GB SSD,IEC60601,Win 11 pro

PC All-In-One, touchscreen, Windows 11 Pro, 23.8", IEC60601, easy to sanitise for medical and dental use.

**Estimated shipping in 10 days**

**SKU:** MOT248PCT

**Categories:** Professional Equipment, All-in-One Medical PCs

### Description

PC All-In-One touchscreen 23.8": computer equipment for viewing medical images must be advanced enough to display faithful graphics and realistic colours. Our touch screen PC with IEC 60601 medical class power supply are designed to meet a wide range of needs, and are especially suitable for displaying top quality images. This device is easy to install and – an indispensable prerequisite for any medical environment – easy to sanitise too. The most important characteristics of this all-in-one medical PC from Tecnomed Italia are its multimedia functionality and its touch screen interface.

#### Technical features PC All-In-One touchscreen 23.8"

Power supply	19 VDC
Flat panel display	23.8? LED IPS touchscreen
Aspect ratio	16:9
Medical power supply	IEC 60601
Resolution	1920 x 1080 Full HD
Processor	Intel Core 5 210H
RAM	16 GB
Total storage capacity	512 GB
Storage media	SSD
Mounting interface	VESA 100
Operating system	Windows 11



Colour	RAL 9003
Touchscreen	Yes
Warranty	24 months
CE	Yes