



PC All-In-One,touchscreen,Windows 11Pro,21,5",IEC60601, mouse, tastiera

PC All-In-One, touchscreen, Windows 11 Pro, 21.5", IEC60601, with mouse & keyboard easy to sanitise for medical and dental use.

Available

SKU: MOT248PCT

Categories: All-in-One Medical PCs

Description

PC All-In-One touchscreen 21.5": computer equipment for viewing medical images must be advanced enough to display faithful graphics and realistic colours. Our touch screen PC with IEC 60601 medical class power supply are de-signed to meet a wide range of needs, and are especially suitable for displaying top quality images. This device is easy to install and – an indispensable prerequisite for any medical environment – easy to sanitise too. The most important characteristics of this all-in-one medical PC from Tecnomed Italia are its multimedia functionality and its touch screen interface.

Technical features PC All-In-One touchscreen 21.5":

Flat screen	21.5" LED Touch screen IPS
Aspect Ratio	16:9
Integrated graphics	Intel UHD Graphics 610
Networking	Wireless LAN – Ethernet LAN – Bluetooth
Power supply	19
Mains adapter	19 VDC /220 VAC



Medical mains adapter	IEC 60601
Resolution	1920 x 1080 Full HD
Microprocessor	Intel Pentium Gold g7400 3.7 Ghz
RAM-HDD	DDR4-SDRAM 8 GB
Total storage	256 GB
Memory support	SSD
Mounting	VESA 100
Input/output	2*USB 3.0, 4*USB 2.0, HDMI
Output	VGA, RJ45, audio
Operating system	Windows 11 professional
Dimensions (w x l x h mm)	546x330x60
Color	RAL 9003
Touchscreen	•
Wireless Keyboard	•
Wireless Mouse	•
Warranty (months)	24
CE	•