



Printed on 31/05/2025

<https://shop.tecnomeditalia.com/en/product/perio-tip-for-scaler-compatible-kavo-piezoled-203/>



Perio tip for scaler, compatible KaVo Piezoled 203

Perio tip for removing subgingival tartar and from deep periodontal pockets. KaVo Piezoled 203 compatible

46,32 €

Available

SKU: PG023K

Categories: Periodontal inserts

Description

Perio tip for scaler PG023K.

Tecnomed Italia inserts for periodontal treatments are extremely thin and able to enter even deep pockets easily, enabling you to achieve excellent and painless results in radicular debridement.

They are also ideal for removing tartar from tight interdental spaces and touching crowns.

The reduced oscillations generated by the low power settings needed for these tips makes treatment with them particularly pleasant for the patient.

Their subgingival spray system eliminates biofilm fragments more effectively.

Perio tip for removing subgingival tartar and from deep periodontal pockets.

Compatibility Perio tip for scaler PG023K:

- KaVo Piezoled 203.