



Perio tip for scaler, compatible KaVo Piezoled 203

Perio tip for removing subgingival tartar and from deep periodontal pockets. KaVo Piezoled 203 compatible

Available

SKU: PG023K

Categories: KaVo

Description

The perio tip for scaler compatible with KaVo Piezoled 203 is a periodontal tip indicated for removing subgingival calculus and calculus from deep periodontal pockets. Its slim shape allows access to deep pockets and difficult-to-reach areas, supporting accurate root debridement during periodontal treatments.

Specifications

- Compatible with KaVo Piezoled 203
- Suitable for use with KaVo Piezoled scaler
- For professional dental use