


Perio tip for scaler, Satelec 1S compatible



Perio tip for removing subgingival tartar and from deep periodontal pockets. Compatible Satelec 1S

Available

SKU: PG099

Categories: Professional Equipment, Satelec - DTE - NSK

Description

Perio tip scaler PG099.

Tecnomed Italia inserts for periodontal treatments are extremely thin and able to enter even deep pockets easily, enabling you to achieve excellent and painless results in radicular debridement.

They are also ideal for removing tartar from tight interdental spaces and touching crowns.

The reduced oscillations generated by the low power settings needed for these tips makes treatment with them particularly pleasant for the patient.

Their subgingival spray system eliminates biofilm fragments more effectively.

Perio tip for removing subgingival tartar and from deep periodontal pockets.

Compatibility Perio tip scaler PG099.

- Satelec 1S.