



Printed on 06/06/2025

<https://shop.tecnomeditalia.com/en/product/perio-tip-for-sine-wave-drive-scaler-compatible-p53l-pl1-p21le/>

Perio tip for SINE WAVE DRIVE scaler, compatible P53L, PL1, P21LE


Titanium perio tip for removing tartar and biofilm. Compatible Woodpecker P53L, EMS PL1, NSK/E P21LE

55,58 €

Available

SKU: PG053LT

Categories: Periodontal inserts



Description

Perio tip for Sine Wave Drive scaler LT.

Tecnomed Italia inserts for periodontal treatments are extremely thin and able to enter even deep pockets easily, enabling you to achieve excellent and painless results in radicular debridement.

They are also ideal for removing tartar from tight interdental spaces and touching crowns.

The reduced oscillations generated by the low power settings needed for these tips makes treatment with them particularly pleasant for the patient.

Their subgingival spray system eliminates biofilm fragments more effectively.

Titanium perio tip for removing tartar and biofilm.

Compatibility Perio tip for Sine Wave Drive scaler LT:

- Woodpecker P53L;
- EMS PL1;
- NSK/E P21LE.