



## Pharmaceutical showcase 1000x580x1970(h)mm, stainless steel, 4 shelves

Available on request

SKU: TM5020SM

Categories: Display Cabinets

### Description

---

The **stainless-steel pharmaceutical cabinet 1000x580x1970 H** is designed to ensure **efficient and hygienic storage of medicines** in **clinics, laboratories and healthcare facilities**.

Made entirely of **AISI 304 stainless steel**, it provides **excellent durability, resistance and ease of cleaning**, guaranteeing long-term performance in demanding medical environments.

The **glass doors** allow quick and clear **visual access to stored medicines**, helping staff check stock levels efficiently.

Inside, the **four adjustable shelves** can be positioned according to different packaging sizes, ensuring a well-organised and functional storage system.

With its sturdy frame and minimalist design, this **pharmaceutical stainless-steel cabinet** combines **practicality, hygiene and reliability**.

It is ideal for **treatment rooms, laboratories and preparation areas**, where organisation and cleanliness are essential.

Smooth surfaces make cleaning operations easy and reduce the accumulation of residues, ensuring compliance with hygiene protocols.

With overall dimensions of **1000x580x1970 H mm**, this **stainless-steel pharmaceutical cabinet** provides a professional, hygienic and space-efficient solution for storing medicines safely.