



Printed on 31/05/2025

<https://shop.tecnomeditalia.com/en/product/piezosurgical-tip-for-scaler-comp-mectron-el1-woodpecker-ul2/>



## Piezosurgical tip for scaler, comp. Mectron EL1, Woodpecker UL2

Piezosurgical tip for separating the Schneider membrane from the edge. Compatible Mectron EL1, Woodpecker UL2

**58,41 €**

**Available**

**SKU:** PG079M

**Categories:** Piezosurgery inserts

### Description

---

Piezosurgical tip for scaler PG079M.

Tecnomed Italia stainless steel piezosurgery inserts are reinforced with a titanium nitride coating to improve reliability and ensure superior accuracy in operations on bone.

Traditional tools like burs and saws can easily provoke permanent damage if they come into contact with nerve fibre.

Tecnomed piezosurgery inserts do not cause damage in the event of occasional or accidental contact with nerve fibre, and because bone is not traumatised by cutting, post-operative recovery is faster.

Compatibility Piezosurgical tip for scaler PG079M.

- Mectron EL1;
- Woodpecker UL2.