



Piezosurgical tip for scaler, comp. Mectron EL1, Woodpecker UL2

Piezosurgical tip for separating the Schneider membrane from the edge. Compatible Mectron EL1, Woodpecker UL2

Available

SKU: PG079M

Categories: Mectron - Woodpecker

Description

Piezosurgical tip for scaler compatible with Mectron EL1 and Woodpecker UL2 is designed for separating the Schneiderian membrane from the margin.

Specifications

- Compatible with Mectron EL1 and Woodpecker UL2
- Suitable for use with compatible Mectron and Woodpecker piezoelectric scalers
- For professional dental use