



Printed on 02/06/2025

<https://shop.tecnomeditalia.com/en/product/piezosurgical-tip-for-scaler-comp-mectron-ex1-woodpecker-uc1/>



## Piezosurgical tip for scaler, comp. Mectron EX1, Woodpecker UC1

Piezosurgical tip for the extraction of teeth in bone inclusions. Compatible Mectron EX1, Woodpecker UC1

**58,41 €**

**Available**

**SKU:** PG082M

**Categories:** Conservative and prosthetic

### Description

---

Piezosurgical tip for scaler PG082M.

The range Conservative & Prosthetics of Tecnomed Italia inserts is ideal for the microinvasive preparation of cavities and for decementation.

Our diamond conservative tips let you target lesions with great accuracy without damaging healthy tissue.

They are also excellent for removing carious dentine without damaging healthy adjacent teeth.

Our inserts for the decementation of implants and bridges are specially designed to help you work more rapidly and effectively.

Compatibility Piezosurgical tip for scaler PG082M.

- Mectron EX1;
- Woodpecker UC1.