



Printed on 31/05/2025

<https://shop.tecnomeditalia.com/en/product/piezosurgical-tip-for-scaler-comp-mectron-op1-woodpecker-us5/>



Piezosurgical tip for scaler, comp. Mectron OP1, Woodpecker US5

Piezosurgical drill for bone remodelling, curettage and bone fragment collection. Compatible Mectron OP1, Woodpecker US5

70,11 €

Available

SKU: PG078M

Categories: Piezosurgery inserts

Description

Piezosurgical tip for scaler PG078M.

Tecnomed Italia stainless steel piezosurgery inserts are reinforced with a titanium nitride coating to improve reliability and ensure superior accuracy in operations on bone.

Traditional tools like burs and saws can easily provoke permanent damage if they come into contact with nerve fibre.

Tecnomed piezosurgery inserts do not cause damage in the event of occasional or accidental contact with nerve fibre, and because bone is not traumatised by cutting, post-operative recovery is faster.

Compatibility Piezosurgical tip for scaler PG078M.

- Mectron OP1;
- Woodpecker US5