



Piezosurgical tip for scaler, comp. Mectron OP1, Woodpecker US5

Piezosurgical drill for bone remodelling, curettage and bone fragment collection. Compatible Mectron OP1, Woodpecker US5

Available

SKU: PG078M

Categories: Mectron - Woodpecker

Description

Piezosurgical tip for scaler, compatible with Mectron OP1 and Woodpecker US5, for bone remodelling, curettage and bone fragment collection in professional dental procedures.

Specifications

- Compatible with Mectron OP1 and Woodpecker US5
- Suitable for use with compatible Mectron and Woodpecker piezoelectric scalers
- For professional dental use