



## Piezosurgical tip for scaler, comp. Mectron OP3, Woodpecker US4

Piezosurgical tip for periodontal osteotomy. Compatible Mectron OP3, Woodpecker US4

**Available**

**SKU:** PG077M

**Categories:** Professional Equipment, Mectron - Woodpecker

### Description

---

Piezosurgical tip for scaler PG077M.

Tecnomed Italia stainless steel piezosurgery inserts are reinforced with a titanium nitride coating to improve reliability and ensure superior accuracy in operations on bone.

Traditional tools like burs and saws can easily provoke permanent damage if they come into contact with nerve fibre.

Tecnomed piezosurgery inserts do not cause damage in the event of occasional or accidental contact with nerve fibre, and because bone is not traumatised by cutting, post-operative recovery is faster.

Compatibility Piezosurgical tip for scaler PG077M.

- Mectron OP3
- Woodpecker US4