



Printed on 31/05/2025

<https://shop.tecnomeditalia.com/en/product/piezosurgical-tip-for-scaler-comp-mectron-ot1-woodpecker-ul4/>



## Piezosurgical tip for scaler, comp. Mectron OT1, Woodpecker UL4

Piezosurgical diamond disc drill for osteotomy in the vicinity of nerves. Compatible Mectron OT1, Woodpecker UL4

**58,41 €**

**Available**

**SKU:** PG081M

**Categories:** Piezosurgery inserts

### Description

---

Piezosurgical tip for scaler PG081M.

Tecnomed Italia stainless steel piezosurgery inserts are reinforced with a titanium nitride coating to improve reliability and ensure superior accuracy in operations on bone.

Traditional tools like burs and saws can easily provoke permanent damage if they come into contact with nerve fibre.

Tecnomed piezosurgery inserts do not cause damage in the event of occasional or accidental contact with nerve fibre, and because bone is not traumatised by cutting, post-operative recovery is faster.

Compatibility piezosurgical tip for scaler PG081M:

- Mectron OT1;
- Woodpecker UL4.