



Printed on 31/05/2025

<https://shop.tecnomeditalia.com/en/product/piezosurgical-tip-for-scaler-comp-mectron-ot2-woodpecker-us3/>



## Piezosurgical tip for scaler, comp. Mectron OT2, Woodpecker US3

Piezosurgical tip for separating the alveolar rim. Compatible  
Mectron OT2, Woodpecker US3

**58,41 €**

**Available**

**SKU:** PG076M

**Categories:** Piezosurgery inserts

### Description

---

Piezosurgical tip for scaler PG076M.

Tecnomed Italia stainless steel piezosurgery inserts are reinforced with a titanium nitride coating to improve reliability and ensure superior accuracy in operations on bone.

Traditional tools like burs and saws can easily provoke permanent damage if they come into contact with nerve fibre.

Tecnomed piezosurgery inserts do not cause damage in the event of occasional or accidental contact with nerve fibre, and because bone is not traumatised by cutting, post-operative recovery is faster.

Piezosurgical tip for separating the alveolar rim.

Compatibility Piezosurgical tip for scaler PG076M:

- Mectron OT2;
- Woodpecker US3.